

Universal Control Arm For use with “Bain” circuits

Add the safety benefit of a Bain non-rebreathing circuit to your existing anesthesia machine.



P/N MA2013
Includes patented safety
pop-off valve

Features:

- Precision manometer to safely measure pressure
- Uses standard 22mm neck rebreathing bags
- Industry standard 22mm circuit port
- Adapts to most anesthesia machines with optional mounting brackets

Includes:

60" long Bain circuit, 1 liter heavy duty rebreathing bag, 3/8" "Z" post, mounting block, screws and wrench for mounting to LEI Medical or Supera stands.

Options:

- Patented safety pressure relief valve (MA2013)
- Mounting brackets for most machines

P/N MA2012 Includes standard pressure relief valve

P/N MA2013 Includes safety pressure relief valve



P/N MA2012
Includes standard pop-off valve

